

Guardian™ Autosampler Probe

The unique design of the robust tip, which combines drip-resistance and built-in particle filtering, helps to prevent cross-contamination during probe movement and blockages in your nebulizer and capillary tubing. Constructed entirely from Ceramic, PEEK and PTFE, the Guardian autosampler probe also provides exceptional resistance to strong acids and solvents.


Features & Benefits

- Robust tip design prevents crushed and damaged tips due to misalignment.
- Drip-resistance prevents cross contamination of samples, especially with oils.
- Built-in particle filtering holds back particulates from blocking the line.
- Completely inert design, Ceramic, PEEK and PTFE construction.
- Interchangeable UniFit™ sample lines available in various IDs (e.g. 0.3, 0.50, 0.75 & 1.0mm)
- Designed to suit Teledyne Cetac®, Agilent®, PerkinElmer®, Shimadzu®, Aim Lab and Thermo Fisher Scientific™ Autosamplers. Thermo PRO is required and sold separately.



Guardian Probe landing page, ordering information, and comparison videos.

[Scan to view](#) >



ASIA PACIFIC

6 Central Boulevard,
Port Melbourne, VIC,
3207, Australia
+61 3 9320 1111
enquiries@geicp.com


AMERICAS

31 Jonathan Bourne
Drive, Unit 7, Pocasset,
MA 02559, USA
800 208 0097
geusa@geicp.com

EUROPE

Friedenbachstrasse
9, 35781 Weilburg,
Germany
+49 6471 3778517
gegmbh@geicp.com

Installed on Leading Autosamplers



P/N 70-803-1787
installed on Cetac



P/N 70-803-2008 installed
on Agilent SPS 4



P/N 70-803-2754 installed
on PerkinElmer S23



P/N 70-803-2849 installed
on PerkinElmer S10

Ordering Information

Guardian Probe with Connecting Line (Red)

P/N	Description
70-803-1957	for Agilent SPS3, SPS4, AIMS, 0.75mm
70-803-1803	for Cetac ASX-200, 500, 800 Series, 0.75mm
70-803-2940	for Cetac ASX-7400, 7600 Series, 0.75mm
70-803-2097	for PerkinElmer S20 Series, 0.75mm
70-803-2851	for PerkinElmer S10, 0.75mm
70-803-2836	for PerkinElmer AS93, 0.75mm
70-803-2872	for Shimadzu AS-10 & AS-20 with Arm Assembly, 0.75mm
70-803-2106	for Thermo iSC-65, 0.75mm

Guardian Probe only

P/N	Description
70-803-2008	for Agilent SPS3, SPS4, AIMS
70-803-1787	for Cetac ASX-200, 500, 800 Series
70-803-2937	for Cetac ASX-7400, 7600 Series
70-803-2029	for Cetac ASX-112FR
70-803-2754	for PerkinElmer S20 Series
70-803-2849	for PerkinElmer S10
70-803-2819	for PerkinElmer AS93
70-803-2860	for Shimadzu AS-10 & AS-20
70-803-2837	for Thermo iSC-65

Probe Connecting Lines

P/N	Description
70-803-1721	Probe Connecting line 1.0mm ID (Green)
70-803-1714	Probe Connecting Line 0.75mm ID (Red)
70-803-1852	Probe Connecting Line 0.5mm ID (Blue)
70-803-1853	Probe Connecting Line 0.3mm ID (Black)

Probe Connecting Lines for ASX-112FR

P/N	Description
70-803-2030	Probe Connecting Line 0.25mm ID (Black)
70-803-2085	Probe Connecting Line 0.18mm ID with EzyFit (Green/Black)

**Note: The above connecting lines are compatible with all Autosamplers (excluding ASX-112FR)*

Visit www.geicp.com for more information, or contact us at enquiries@geicp.com