

Glass DC Nebulizers

SeaSpray

| | |
|------------------|--------------|
| 2mL/min uptake | A13-07-USS2 |
| 400µL/min uptake | A13-07-USS04 |

Conikal

| | |
|----------------|------------|
| 2mL/min uptake | A13-07-UC2 |
| 1mL/min uptake | A13-07-UC1 |

Slurry

| | |
|----------------|------------|
| 4mL/min uptake | A13-07-US6 |
|----------------|------------|



DC Nebulizers for HF

DuraMist

| | |
|----------------|------------|
| 1mL/min uptake | A13-07-DM1 |
|----------------|------------|

OpalMist

| | |
|----------------|-------------|
| 2mL/min uptake | A13-07-PFA2 |
|----------------|-------------|

VeeSpray

| | |
|----------------------|------------|
| 0.6 - 3mL/min uptake | A13-07-CV6 |
|----------------------|------------|



IsoMist XR

Programmable Temperature Cyclonic Spray Chamber

| | |
|----------------|------------|
| IsoMist XR Kit | KT-1137-XR |
|----------------|------------|

Comprises:

| | |
|--|-------------|
| Glass Twister Spray Chamber | 21-809-3186 |
| Assembly (ROL) | 20-809-5176 |
| PuraFlex Transfer Tube | 21-809-4038 |
| Metal Ball Joint Clip | JC-13 |
| UniFit Right Angle Connector 2.0mm (PKT 3) | GAZ-04-2-90 |




BC Spray Chambers

Tracey BC

| | |
|--|-------------|
| 30mL Cyclonic | 20-809-5316 |
| * Direct replacement for 20-808-8882HE | |

Tracey BC PEEK

| | |
|--------------------------------------|-------------|
| 30mL Cyclonic | 20-809-4795 |
| * Direct replacement for 20-809-2507 | |
| * No internal resurfacing required | |



NEW*



NEW*

Spray Chambers

Twister

| | |
|---------------|---------------|
| 50mL Cyclonic | 20-809-9199HE |
|---------------|---------------|

Tracey

| | |
|---------------|---------------|
| 50mL Cyclonic | 20-808-8882HE |
|---------------|---------------|



HydraMist Simultaneous Cold Vapor/ Pneumatic Nebulization Spray Chamber

| | |
|-------------|---------|
| Hydride Kit | KT-1157 |
|-------------|---------|



Spray Chambers for HF

Twister TFE Spray Chamber

| | |
|---------------|-------------|
| 50mL Cyclonic | 20-809-3452 |
|---------------|-------------|

Tracey TFE Spray Chamber

| | |
|---------------|-------------|
| 50mL Cyclonic | 20-809-2507 |
|---------------|-------------|




Torches

D-Torch

SVDV/VDV 30-808-3560 


Comprises:

Base and Inner Tube 31-808-3568
 Quartz Outer Tube 31-808-3570
 Retaining Ring 31-808-3566
 Injector Ferrule 31-808-2962
 Injector Adaptor 31-808-3578 



D-Torch

RV 30-808-3590

Comprises:

Base and Inner Tube 31-808-3568 
 Quartz Outer Tube 31-808-3582
 Retaining Ring 31-808-3566
 Injector Ferrule 31-808-2962
 Injector Adaptor 31-808-3578

D-Torch Options:

Ceramic Outer Tube for SVDV/VDV 31-808-3580 
 Ceramic Outer Tube for RV 31-808-3586
 Quartz Outer Tube for SVDV/VDV 31-808-3570 
 Quartz Outer Tube for RV 31-808-3582
 O-ring Kit 70-803-1500
 Ignition Enhancer 31-808-4020

Quartz Torch

SVDV/VDV 31-808-3557 

RV 31-808-3556

Retaining Ring 31-808-3563 

Pre-Optic Window Axial 70-803-1968

Pre-Optic Window Radial 70-803-1969

O-ring for Window Axial (PKT 5) 70-803-2794 

Pre-optics cone 70-803-2751

Torches

(continued)

D-Torch Injector Options

Quartz Injector 0.8mm 31-808-3738 

Quartz Injector 1.4mm 31-808-3583

Quartz Injector 1.8mm 31-808-3576

Quartz Injector 2.4mm 31-808-3584

Alumina Injector 1.5mm 31-808-3725 

Alumina Injector 1.8mm 31-808-3585

Alumina Injector 2.4mm 31-808-3722

Injectors for Agilent FDT

Quartz Injector 0.8mm 31-808-4357

Quartz Injector 1.4mm 31-808-4354 

Quartz Injector 1.8mm 31-808-4355

Quartz Injector 2.4mm 31-808-4356

Alumina Injector 0.8mm 31-808-4530

Alumina Injector 1.4mm 31-808-4532 

Alumina Injector 1.8mm 31-808-4358

Alumina Injector 2.4mm 31-808-4534

RF Coils

Silver 70-900-5100S 

Gold 70-900-5100G

Copper 70-900-5100C

RF Coil Installation Kit 70-900-5100T

RotaValve

RotaValve 70-803-1738 

Valve Tubing Kit 70-803-1869 

Sample Loop 1.0mL 70-803-1860 

Nebulizer Connector UniFit - 1/4-28 Flangless 70-803-1888 

Common Accessories

Elegra Argon Humidifier

Elegra Argon Humidifier 70-803-1268



Guardian In-Line Filters

Guardian DC in-line non-return gas filter 70-803-1942

Guardian in-line sample filter 70-803-1108


Eluo Adaptor for Inline Filter 70-803-1160



Trident Internal Standard Kits

Trident CT Standard Kit 60-703-1179

Standard Kit for HF 60-808-1150



Magnifier Inspection Tool 70-803-1923



Autosampler Probes


(continued)

For SPS 3/SPS 4

Guardian Probe 0.75mm
Probe Connecting Line (Red) 70-803-1957

PTFE Encapsulated Carbon Fibre Probe
0.75mm ID 70-803-0908

PTFE Encapsulated Carbon Fibre Probe
1.0mm ID 70-803-0853



Autosampler Probes (See page 22-23 for complete list of probes)

For Cetac ASX-200/500/800 Series

Guardian Probe 0.75mm
Probe Connecting Line (Red) 70-803-1803

Polyimide sheathed Autosampler Probe
0.5mm ID 60-808-1186L

PTFE Encapsulated Carbon Fibre Probe
0.5mm ID 70-803-0784

PTFE Encapsulated Carbon Fibre Probe
0.75mm ID 70-803-1443

PTFE Encapsulated Carbon Fibre Probe
1.0mm ID 70-803-0793



For Cetac ASX-7400, 7600 Series

Guardian Probe 0.75mm
Probe Connecting Line (Red) 70-803-2940

