

Glass DC Nebulizers

SeaSpray

400µL/min uptake A23-1-USS04

MicroMist

400µL/min uptake A23-1-UM04

Conikal

1mL/min uptake A23-1-UC1



Quartz DC Nebulizers

SeaSpray

1mL/min uptake A23-1-QSS1



DC Nebulizers for HF

DuraMist

400µL/min uptake A23-1-DM04



OpalMist

400µL/min uptake A23-1-PFA04



Spray Chambers

Twister

50mL Cyclonic 20-809-4417HE

50mL Cyclonic 20-809-4560

With Matrix Port



Tracey

50mL Cyclonic 20-809-4419HE



Spray Chambers


Cinnabar

20mL Cyclonic 20-809-4480HE



High Efficiency Sample Introduction System

HE-SIS for PerkinElmer NexION KT-1184



Torches

D-Torch

30-808-3927

Comprises:

- Base and Inner Tube 31-808-3928
- Quartz Outer Tube 31-808-3942
- Retaining Ring 31-808-2818
- Injector Adaptor 31-808-3980
- Injector Ferrule 31-808-2962
- Fastener Assembly 31-808-4130
- Viton O-rings (PKT 10) 70-V-011
- PTFE Coated Viton O-rings (PKT 10) 70-VT-012



Quartz Injectors

Tapered 1.0mm 31-808-3423

Tapered 1.8mm 31-808-4005

Tapered 2.0mm 31-808-4006

Tapered 2.5mm 31-808-3302



HF Resistant Injectors

Tapered Alumina 2.4mm 31-808-3205

Sapphire 2.4mm 31-808-3378



Cones

Sampler Cones

Nickel Sampler PE3011-Ni



Platinum Sampler PE3013-Pt

Platinum Sampler (18mm) PE3013A-Pt

Platinum Sampler (12mm) PE3013C-Pt

Platinum Sampler PE3013-Pt

Boron Free PE3013-Pt-BF



Skimmer Cones

Nickel Skimmer PE4212-Ni

Platinum Skimmer PE4214-Pt

Nickel Hyper Skimmer PE4015-Ni

All Nickel Hyper Skimmer PE4016-Ni



O-ring for NexION Hyper Skimmer Cone (pack of 5)

PE5125



ConeGuard Thread Protectors

Sampler Cone Guard 70-803-1024
For NexION



Gasket for Sampler Cone (pack of 5)

PE5111



Common Accessories

Elegra Argon Humidifier

Elegra Dual Argon Humidifier 70-803-1473



Guardian In-Line Filters

Guardian DC in-line non-return gas filter 70-803-1942

Guardian in-line sample filter 70-803-1108


Eluo Adaptor for Inline Filter 70-803-1160



Trident Internal Standard Kits

Trident CT Standard Kit 60-703-1179

Standard Kit for HF 60-808-1150



Magnifier Inspection Tool 70-803-1923



Autosampler Probes (See page 22-23 for complete list of probes)

For PerkinElmer S20 Series

Guardian Probe 0.75mm
Probe Connecting Line (Red) 70-803-2097



For Cetac ASX-200/500/800 Series

Guardian Probe 0.75mm
Probe Connecting Line (Red) 70-803-1803

Polyimide sheathed Autosampler
Probe 0.5mm ID 60-808-1186L

PTFE Encapsulated Carbon Fibre
Probe 0.5mm ID 70-803-0784

PTFE Encapsulated Carbon Fibre
Probe 0.75mm ID 70-803-1443

PTFE Encapsulated Carbon Fibre
Probe 1.0mm ID 70-803-0793

