

Glass DC Nebulizers

SeaSpray

2mL/min uptake A21-07-USS2

MicroMist

600µL/min uptake A21-07-UM06

Slurry

4mL/min uptake A21-07-US6

Conikal

2mL/min uptake A21-07-UC2



Quartz DC Nebulizers

SeaSpray

1mL/min uptake A21-07-QSS1



DC Nebulizers for HF

DuraMist

1mL/min uptake A21-07-DM1



OpalMist

2mL/min uptake A21-07-PFA2



IsoMist XR

Programmable Temperature Cyclonic Spray Chamber

IsoMist XR Kit KT-1012-XR

Comprises:

Glass Twister Spray Chamber 21-809-3185

Assembly (LOL) 20-809-5175

IsoMist Mount 21-809-2773

Adaptor GAZ08-12 90deg 21-809-2871

Clamping Screw 21-808-0953

UniFit Right Angle Connector (PKT 3) GAZ-04-90



Spray Chambers

Tracey

50mL Cyclonic 20-809-0434HE



Twister

50mL Cyclonic (Aqueous) 20-809-0208HE
With 7mm baffle



Twister

50mL Cyclonic (Organics) 20-809-0446HE
With 4mm baffle



Twister

50mL Cyclonic (Aux. Port) 20-809-0004HE



Cinnabar

20mL Cyclonic 20-809-9121HE



HydraMist Simultaneous Cold Vapor/ Pneumatic Nebulization Spray Chamber

HydraMist Kit KT-1162



Spray Chambers for HF

Tracey TFE

50mL Cyclonic 20-809-2510

Plugs directly into cassette mount




Torches

D-Torch 30-808-3300 

Comprises:

- Base and Inner Tube 31-808-3292
- Quartz Outer Tube 31-808-3294
- Retaining Ring 31-808-2818

D-Torch Options

- Ceramic Outer Tube 31-808-2990 
- Copper Foils (PKT 12) for Cermic Outer Tube 31-808-3356
- O-Ring Kit 70-803-0762

Quartz Torch

- Without slots 31-808-3254 
- With single slot 31-808-3255
- With 3 slots 31-808-3256

Quartz Torch Options:

- Copper Foil for Quartz Torch or D-Torch Quartz Outer Tube (pack of 12) 31-808-3289


Tapered Quartz Injectors

- Quartz Injector 0.8mm 31-808-3257 
- Quartz Injector 2.0mm 31-808-3258
- Alumina Injector 1.2mm 31-808-3287
- Alumina Injector 2.0mm 31-808-3286 

Sapphire Injectors

- Sapphire Injector 2mm 31-808-3528

Purge Windows

- Axial 17mm disk 31-808-0551 
- Axial 32mm disk 31-808-2649
- Radial (Post Nov. 2004) 31-808-2788

Common Accessories


Elegra Argon Humidifier

- Elegra Argon Humidifier 70-803-1269 



Common Accessories

(continued)

Guardian In-Line Filters

- Guardian DC in-line non-return gas filter 70-803-1942 
- Guardian in-line sample filter 70-803-1108
- Eluo Adaptor for Inline Filter 70-803-1160



Trident Internal Standard Kits

- Trident CT Standard Kit 60-703-1179 
- Standard Kit for HF solutions 60-808-1150 



Magnifier Inspection Tool 70-803-1923 

Autosampler Probes (See page 22-23 for complete list of probes)

Suitable for PE S10 or AS93+

- Guardian Probe 0.75mm Probe Connecting Line (Red) (*for PE S10 only) 70-803-2851 
- Guardian Probe 0.75mm Probe Connecting Line (Red) (*for AS93+ only) 70-803-2836 
- PTFE Encapsulated Carbon Fibre Probe 0.25mm ID 70-803-1071 
- PTFE Encapsulated Carbon Fibre Probe 0.5mm ID 70-803-1440
- PTFE Encapsulated Carbon Fibre Probe 1.0mm ID 70-803-0816

Suitable for Cetac ASX-200/500/800 Series

- Guardian Probe 0.75mm Probe Connecting Line (Red) 70-803-1803 
- Polyimide sheathed Autosampler Probe 0.5mm ID 60-808-1186L
- PTFE Encapsulated Carbon Fibre Probe 0.5mm ID 70-803-0784 
- PTFE Encapsulated Carbon Fibre Probe 0.75mm ID 70-803-1443
- PTFE Encapsulated Carbon Fibre Probe 1.0mm ID 70-803-0793

For Cetac ASX-7400, 7600 Series

- Guardian Probe 0.75mm Probe Connecting Line (Red) 70-803-2940 