

Glass DC Nebulizers

SeaSpray

2mL/min uptake A31-07-USS2

MicroMist

600µL/min uptake A31-07-UM06

Conikal

2mL/min uptake A31-07-UC2

Slurry

4mL/min uptake A31-07-US6



DC Nebulizers for HF

DuraMist

1mL/min uptake A31-07-DM1

OpalMist

2mL/min uptake A31-07-PFA2



IsoMist XR

Programmable Temperature Cyclonic Spray Chamber

IsoMist XR Kit KT-1166-XR

Comprises:

Glass Twister Spray Chamber 21-809-3185

Assembly (LOL) 20-809-5175

PuraFlex Transfer Tube 21-809-4426

UniFit Right Angle Connector (PKT 3) GAZ-04-90



Spray Chambers

Tracey

50mL Cyclonic 20-809-2665HE

Twister

50mL Cyclonic 20-809-2951HE

Straight Ball Joint Adaptor

21-809-2723

Sheath Gas Adaptor

31-808-3891

HydraMist Simultaneous Cold Vapor/ Pneumatic Nebulization Spray Chamber

Hydride Kit KT-1160



Spray Chambers for HF

Tracey TFE

50mL Cyclonic 20-809-9198




Torches

E-Torch 30-808-5163

Comprises:

- Body 31-808-4425
- Quartz Tube Set 31-808-5164
- Injector Adaptor 31-808-4582
- Retaining Ring 31-808-4498




D-Torch Options:

- Ceramic Tube Set 31-808-5165
- FKM O-Rings (PKT 10) 31-808-4381
- Silicone O-rings 22X1.8S70 (PKT 10) 70-803-6475



- Viton O-rings for FDT/SDT ball joints (PKT 10) 70-V-011

- Injector Adaptor S19 31-808-5181



Tapered Injectors for E-Torch

- Quartz Injector 0.75mm 31-808-3511
- Quartz Injector 1.0mm 31-808-2833
- Quartz Injector 1.5mm 31-808-2834



Capillary Injectors for E-Torch

- Quartz Injector 1.0mm 31-808-3510
- Quartz Injector 2.0mm 31-808-2836
- Quartz Injector 2.5mm 31-808-3246
- Quartz Injector 3.0mm 31-808-3247
- Alumina Injector 1.0mm 31-808-3279
- Alumina Injector 2.0mm 31-808-2835
- Sapphire Injector 2.0mm 31-808-2977



Torches

(continued)

D-Torch 30-808-2877

Comprises:

- Base and Inner Tube 31-808-2932
- Quartz Outer Tube 31-808-2878



D-Torch Options:

- Ceramic Outer Tube 31-808-3836
- PTFE Coated
- Viton O-rings (PKT 10) VT-007
- Viton O-rings (PKT 10) V-026



- O-ring kit for Thermo 6000/7000/PRO D-Torch 70-803-1921



- O-ring and Seal kit for Thermo Torch & Windows 70-803-1550

- Retaining Ring 31-808-2707



- Injector Adaptor 31-808-2895

- Injector Ferrule 31-808-2897

Quartz Torch

- EMT Design 31-808-2810



Tapered Injectors for D-Torch & EMT Torch

- Quartz Injector 0.75mm 31-808-3511
- Quartz Injector 1.0mm 31-808-2833
- Quartz Injector 1.5mm 31-808-2834



Capillary Injectors for D-Torch & EMT Torch

- Quartz Injector 1.0mm 31-808-3510
- Quartz Injector 2.0mm 31-808-2836
- Quartz Injector 2.5mm 31-808-3246
- Quartz Injector 3.0mm 31-808-3247
- Alumina Injector 1.0mm 31-808-3279
- Alumina Injector 2.0mm 31-808-2835
- Sapphire Injector 2.0mm 31-808-2977



RF Coils

RF Coil for PRO

Silver 70-900-4006S



Gold 70-900-4006G



RF Coil Installation Kit 70-900-4006T



Common Accessories

Elegra Argon Humidifier

Elegra Argon Humidifier 70-803-1268



Guardian In-Line Filters

Guardian DC in-line non-return gas filter 70-803-1942

Guardian in-line sample filter 70-803-1108


Eluo Adaptor for Inline Filter 70-803-1160



Trident Internal Standard Kits

Trident CT Standard Kit 60-703-1179

Standard Kit for HF 60-808-1150



Magnifier Inspection Tool 70-803-1923



Autosampler Probes (See page 22-23 for complete list of probes)

For Cetac ASX-200/500/800 Series

Guardian Probe 0.75mm
Probe Connecting Line (Red) 70-803-1803

Polyimide sheathed Autosampler
Probe 0.5mm ID 60-808-1186L

PTFE Encapsulated Carbon
Fibre Probe 0.5mm ID 70-803-0784

PTFE Encapsulated Carbon
Fibre Probe 0.75mm ID 70-803-1443

PTFE Encapsulated Carbon
Fibre Probe 1.0mm ID 70-803-0793



For Cetac ASX-7400, 7600 Series

Guardian Probe 0.75mm
Probe Connecting Line (Red) 70-803-2940

